CareFusion

# BabyCO Meter

Operating Manual – English

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## Contents

Overview .....	9
<b>I. Introduction - English.....</b>	<b>9</b>
Cigarette Consumption and CO Measurements .....	10
Package Contents .....	12
Contraindications .....	13
Warnings and Cautions .....	13
Intended Use .....	14
Operation .....	15
PC Connection .....	17
Calibration .....	18
Battery Life .....	22
Battery Replacement.....	22
Internal Battery Expiry .....	23
Sensor Expiry .....	24
Cleaning .....	25
Servicing .....	26
Trouble Shooting Information.....	26
Safety Designation per IEC 60601-1 .....	27
Electromagnetic Compatibility (EMC) to EN60601-1-2:2007 .....	28
Symbols .....	33
Specifications .....	34
Consumables / Accessories .....	35
Customer contact information .....	37
Übersicht .....	41
<b>II. Einleitung - Deutsch .....</b>	<b>41</b>
Zigarettenkonsum und CO-Messungen .....	42

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## Operating manual

### Overview

The BabyCO meter is a hand held battery operated device used to measure the concentration of carbon monoxide, CO, on the breath and calculates the percentage of carboxyhaemoglobin, FCOHb, in the foetal blood.

It is accurate, easy to use, and has many features designed to simplify its operation.

These include:	Auto zero function	Smoking level colour
light indicators		Simple calibration Serial interface to PC

**Please note: It is recommended that this unit be calibrated upon receipt.**

### I. Introduction - English

The BabyCO meter is based on an electrochemical fuel cell sensor, which works through the reaction of carbon monoxide (CO) with an electrolyte at one electrode, and oxygen (from ambient air) at the other. This reaction generates an electrical current proportional to CO concentration. Output from the sensor is monitored by a microprocessor, which detects peak expired concentrations of alveolar gas. This is then converted to % foetal carboxyhaemoglobin (%FCOHb) using the relationship  $\%FCOHb = 0.313 * Expired\ CO$

(ppm). This relationship is derived from data supplied by the Centre Hospitalier d'Arras, Maternité Georges Pernin, Arras, France<sup>1</sup>. The results are displayed on a clear LCD display.

Warning lights are provided to give an instant indication of the smoking level.

The warning light levels are user adjustable when the unit is connected to a PC running COBRA software.

## References

1. Gomez C, Berlin I, Marquis P, Delcroix M.

**Expired air carbon monoxide concentration in mothers and their spouses above 5 ppm is associated with decreased foetal growth.**

Preventative Medicine 40 (2005) 10-15

## Cigarette Consumption and CO Measurements

The BabyCO meter provides a simple screening test for cigarette consumption for use in anti-smoking clinics and all smoking cessation programs.

Measurement of carboxyhaemoglobin has been well validated as an indirect measure of cigarette consumption and is widely used in smoking cessation programs.

Typical values for carboxyhaemoglobin and expired CO in smokers together with the alarm light status, are given below:

CO(ppm)	Cigarette consumption	Indicator
0 – 6	Non smoker	Green
7 – 10	Light smoker	Amber
11 - 20	Heavy smoker	Red

20+	Heavy smoker	Flashing red with alarm
-----	--------------	-------------------------

Please note that some urban areas may have high environmental levels of CO. This can cause a rise in exhaled CO of a few ppm above that which is normally present on the breath. In these cases it is possible for a non-smoker to appear at the bottom of the 'light smoker' range (7 – 10 ppm).

The above table applies when the unit is measuring expired CO with the switch in middle position. When measuring FCOHb, with the switch in the top position, the red indicator will illuminate whenever the expired CO measurement is 7 ppm and above.


## References

1. **Carbon monoxide in breath in relation to smoking and carboxyhemoglobin levels.**  
Wald NJ, Idle M, Boreham J, Baily A  
Thorax 1981; 36, 366-369
2. **Definition of a reliable threshold value for detecting current smokers by CO measurement**  
Marino Luigi; Latini Roberto; Barbano Gina; Bazzerla Giorgio;  
De Luca Anita, Nardini Stefano - Respiratory and TB Unit-  
General Hospital- Via forlanini, 71-I-31029-Vittorio Veneto (TV-  
ITALY).

## Package Contents

The BabyCO meter is supplied with a carrying case containing the following items:

1. BabyCO meter (Cat No. 36-BC01-STK)
2. 22 mm mouthpiece adapter (Box of 10 Cat No.36-PSA2000)
3. 22 mm reducing connector for calibration (Cat No.36-MEC1007)
4. PP3 Battery (Cat No.36-BAT1002)
5. Disposable mouthpieces (Box 250 Cat No.36-PSA1200)
6. Calibration tool (Cat No. 36-MEC1184)

## Contraindications

There are no known contraindications, however it should be noted:

The BabyCO should not be used on individuals that may cause discomfort exhaling at a sustained speed or exacerbate any discomfort at the mouth by exhaling into a mouthpiece

Pregnant individuals should hold their breath for as long as they feel comfortable to do so before exhaling into the device

## Warnings and Cautions

The following terms are used as follows in this manual

**Caution:** Possibility of injury or serious damage

**Warning:** conditions or practices that could result in personal injury

**Please Note:** Important information for avoiding damage to the instrument or facilitating operation of the instrument.

**Note:** The device should only be used by trained and qualified personnel.

**CAUTION:** Read the manual before use

**CAUTION:** Do not attempt to charge the batteries, connect improperly, or dispose of in a fire as there is possibility of leakage or explosion. Follow manufacturer's recommendation for proper disposal.

**WARNING:** The instrument is not suitable for use in the presence of explosive or flammable gases, flammable anaesthetic mixtures or in oxygen rich environments.

**CAUTION:** Mouthpieces are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may degrade materials and lead to an incorrect measurement.

**PLEASE NOTE:** The product you have purchased should not be disposed of as unsorted waste. Please utilise your local WEEE collection facilities for the disposal of this product.

**PLEASE NOTE:** Degree of protection against Ingress of Water is IPX0.

**CAUTION:** When you connect the BabyCO to other equipment, always make sure that the whole combination complies with the international safety standard IEC 60601-1 for medical electrical systems. During measurements, connect the BabyCO only to computers that comply with IEC/EN 60601-1 / UL 60601-1.

## Intended Use

The BabyCO is intended to measure the amount of Carbon Monoxide (CO) that can be expired from the breath. The expired CO is measured in parts per million and this can be used to calculate using an algorithm the percentage of carboxyhaemoglobin, (%FCOHb), in the foetal blood.


The BabyCO is used in smoking cessation clinics, GP's surgeries and antenatal clinics for instructional purposes and to check on the client's progress and compliance.

The BabyCO is designed for use by clinicians and health care professionals

## Operation

For accurate results the BabyCO meter should be used at room temperature. If the instrument has been stored in cool or hot conditions then allow time to reach room temperature prior to use.


Install the PP3 battery by sliding open the battery cover, clipping the battery in place and replacing the cover. Insert the mouthpiece adapter into the BabyCO meter and then insert a disposable mouthpiece into the adapter. Turn the unit on by pushing the slide switch to the middle position (CO). The display will show the software version number momentarily whilst correction for ambient levels is executed. During this time the unit must not be exposed to any source of CO. The buzzer will then sound and the display will change to:

Upon hearing the buzzer instruct the subject to inspire fully and hold their breath for as long as is comfortable. The breath-holding period is recommended to allow time for equilibration of alveolar gas.

The subject should then seal their lips around the mouthpiece and exhale slowly and fully.

After emptying the lungs completely the display will show the expired CO concentration in the expired air in parts per million (ppm):

Provided the subject is not too distressed after seeing the breath CO measurement, then the calculated foetal carboxyhemoglobin level may be displayed by pushing the slide switch to the top position (FCOHb):

The unit must not be used, however, for one second after the unit is first switched on i.e. before the display shows 'BLOW'. Expired alveolar gas is then entrapped between sensor and mouthpiece valve. The display reading will rise to a plateau over the course of several seconds. The final value will be held until the unit is turned off and represents parts per million CO or %FCOHb depending on the slide switch position.

The lights will come on according to the table on page 10.


**Important note:**

**Before repeating a measurement the unit must be turned off, and the mouthpiece and adapter removed for at least 1 minute. This is to allow re-equilibration with ambient air and to dry the surface of the sensor. Visually inspect that all**

**moisture has evaporated from the surface of the sensor before reuse.**

If the unit is switched on again too quickly after use there may be a response to residual expired carbon monoxide from the previous test.

In this case the display will show:

If this is displayed then turn the unit off, remove the mouthpiece adapter, and expose to ambient air for 2 minutes before repeating the test.

**Note:** If this warning appears again after following the above procedure then turn off and leave the sensor exposed to ambient air for a further 3 minutes. If the same message appears again then this indicates possible contamination of the sensor with a solvent.

In this case remove all sources of solvent from the vicinity of the sensor and expose to ambient air for 24 hours before switching on again.

## **PC Connection**

The BabyCO meter may be connected to the serial port of a PC running COBRA (CO Breath Analysis) software.

This software allows the measurements to be read by a PC and automatically entered onto a pre-defined report for subsequent printing and filing. It also allows the CO level indicators and the

breath holding countdown timer to be configured to individual requirements.

**Note:** The BabyCO should only be connected to a computer that is manufactured in accordance with EN 60601-1.


**Note:** Keep the PC out of reach of the patient at all times.

## Calibration

Calibration will remain stable to within 2% over one month and typically to within 10% over 6 months. CareFusion supplies calibration gas (20ppm CO in air) and recommends that the unit be recalibrated on a 6 monthly basis. See page 35 for calibration accessories.


The gas is supplied in convenient, disposable, aluminium aerosol cans containing 20 litres of gas (Cat. No.36-MCG020).

To carry out the calibration locate the calibration switch on the right-hand side of the instrument as shown below.

Push the slide switch to the CO - PPM position and wait for the unit to display zero.


Screw the control valve firmly onto the cylinder and connect the gas supply as shown below:

The plastic tubing supplied with the gas should be pushed firmly over the reducing connector.


Slowly turn the control knob anti-clockwise until the ball in the flow indicator is between the two marks. This will then supply a gas flow of approximately 0.25 l/min. Apply this flow for 25 seconds and if the meter does not read 20ppm use the calibration tool to push the calibration switch.

The unit will then beep 3 times, store the new calibration value, and display the following:

The gas supply should then be turned off.

If the signal from the CO sensor is too low a new calibration value will not be stored and the display will show:


The most likely cause for this is an expired cell but may also be caused by depressing the calibration switch with no gas applied. Ensure that the concentration of calibration gas is correct (20ppm), that the connections to the gas cylinder are secure, and that the gas cylinder is not empty, and then repeat the calibration procedure.

If the above message is repeated, the BabyCO meter must be returned to CareFusion Ltd, or an authorised agent, for sensor replacement.

Sensor life ranges from 2 to 5 years and depends upon both the amount of exposure to CO and other gases, particularly solvents such as alcohol and cleaning fluids.

If the signal from the CO sensor is too high, a new calibration value will not be stored and the following will be displayed when calibration is attempted:

The most likely cause for this is using an incorrect concentration of calibration gas. Ensure that the concentration of calibration gas is correct (20ppm), that the connections to the gas cylinder are secure, and then repeat the calibration procedure.

In order to stop any gas leak from the can after calibration, ensure that the control knob is tightened firmly.

### **Important Notes:**


- Only certified calibration gas from a reputable source should be used.
- Ensure that no CO is present on the sensor for 3 minutes before starting the calibration procedure.
- Ensure that the instrument and gas cylinder have stabilised at room temperature before calibrating.

- In order to stop any gas leak from the can after calibration, ensure that the control knob is tightened firmly.

## Battery Life

Battery life is approximately 30 hours of continuous use.

When the battery has approximately 1 hour of useful life left the alarm will sound momentarily after the unit is first switched on and the following message will be displayed:

When the battery has completely expired the above will be displayed continuously and the battery must be replaced.

## Battery Replacement

Locate the sliding cover situated on the rear of the unit, towards the bottom of the device.

Place your thumb over the round thumb indent, press gently and slide the cover to the right to remove it from the unit.

Lift the old battery out, and holding the battery terminal by the plastic body, pull it off the old battery.

Plug the new battery into the battery terminal, taking care that the correct polarity is observed.

Push the battery back into the battery holder, and replace the battery cover onto the guides. Slide the battery cover to the left until it is fully home.

**Note:** Please remove the battery if the meter is likely to be unused for some time.

**CAUTION:** Do not open the battery cover when the device is turned on.


**CAUTION:** The operator should not touch the contacts of the battery and the patient at the same time.

**Please Note:** Dispose of the waste battery in accordance with EU Waste Battery Regulations.

### **Internal Battery Expiry**

The BabyCO meter has an internal battery with a life of approximately 10 years. This battery supplies the sensor signal conditioning circuit continuously to ensure instant start up.

When the battery has expired, the following warning message will be displayed:

and the alarm will sound when the unit is first turned on.

When this occurs, the BabyCO meter must be returned to CareFusion, or an authorised agent, for battery replacement.

## **Sensor Expiry**

Sensor life ranges from 2 to 5 years and depends upon both the amount of exposure to CO and other gases and solvents such as alcohol and cleaning fluids.

When the sensor has expired it will become impossible to obtain a correct calibration. When this occurs, the Baby CO must be returned to CareFusion or an authorised agent, for sensor replacement.

## Cleaning

Disinfection of contaminated parts is only effective after having them carefully pre-cleaned. CareFusion recommend the tested solution of PeraSafe sterilizing powder (36-SSC5000A) for pre-cleaning and disinfection. If a different solution is used, please follow the given manufacturer's instructions.

**CAUTION:** Switch off the device and always unplug your BabyCO meter from the computer before cleaning.

When using the disposable cardboard mouthpiece with one-way valve under the prerequisite that the patient was instructed only to exhale, the mouthpiece adapter and the CO meter surface have to be cleaned.

The mouthpiece adapter may be cleaned and sterilised by immersion in Perasafe (36-SSC5000A) or other cold sterilising solutions.

Exposed surfaces of the meter may be wiped with a clean, damp cloth or a cloth that has been immersed in a sterilising solution. It is recommended that this procedure is performed after each use and all used cardboard mouthpieces discarded.

**Important note:** The sensor surface must not be wiped with any aqueous solutions and **must not** be exposed to solvents e.g. alcohol or permanent damage may result

**CAUTION:** Do not attempt to wash or immerse the BabyCO meter in water or cleaning fluid, as there are electronic components inside that will be permanently damaged.

**Important note:** Cardboard mouthpieces must be disposed off immediately after the use.

If there are changes on the material surfaces of either the unit or mouthpiece adaptor (cracks, brittleness) the respective parts must be disposed off.

## Servicing

If your unit requires service or repair please see page 37 for contact details. A full service manual including circuit diagram and parts list is available upon request.

## Trouble Shooting Information

Should you encounter problems operating your BabyCO meter, please consult the table below:

Problem	Possible Cause	Solution
BabyCO cannot be switched on or "bat" is displayed	Batteries are flat	Change the battery
Battery life is shorter than expected	Unit not being switched off	Turn the unit off after use
"CEL" is displayed	Fuel cell is depleted	Return unit for cell replacement
"Err" is displayed	Gas cylinder empty	Check valve on cylinder for contents and replace cylinder if necessary Check the calibration value on the cylinder is 20ppm.
"bt2" is displayed	Internal battery has expired	Return unit for battery replacement
"gAS" is displayed	Fuel cell is depleted	Return unit for cell replacement
	Fuel cell is contaminated by residual gas or solvents	Allow unit to be exposed to clean air. Turn unit ON and await message to disappear, if not return unit for cell replacement.

## Safety Designation per IEC 60601-1

<b>Type of protection against electrical shock</b>	Internally powered Equipment
<b>Degree of protection against electrical shock</b>	Type B applied part
<b>Power Equipment</b>	Battery type: PP3 battery, Alkaline Manganese Dioxide, 9.0V, 550mAh
<b>Battery life:</b>	2000 tests.
<b>Degree of Electrical connection between equipment and Patient</b>	Equipment designed as non-electrical connection to the patient.
<b>Degree of mobility</b>	Transportable
<b>Mode of operation</b>	Continuous

## Classification according to IEC 60601-1

BabyCO Meter

Applied part, type B

**WARNING:** No modification of this equipment is allowed.

**WARNING:** Do not connect devices that are not specified as part of the system.

**NOTE:** When you connect other equipment to the unit, always make sure that the whole combination complies with the international safety standard IEC 60601-1-1 for medical electrical systems. When connecting to a PC with the supplied serial cable then the PC must be EN 60601-1 / UL 60601-1 compliant.

**WARNING:** The user must not touch any voltage-carrying parts and the patient at the same time.

## **Electromagnetic Compatibility (EMC) to EN60601-1-2:2007**

**WARNING:** use of portable phones or other radio frequency (RF) emitting equipment near the system may cause unexpected or adverse operation

The BabyCO has been tested to EN60601-1-2:2007, regarding it's ability to operate in an environment containing other electrical/electronic equipment (including other medical devices).

The purpose of this testing is to ensure that the BabyCO is not likely to adversely affect the normal operation of other such equipment and that other such equipment is not likely to adversely affect the normal operation of the BabyCO.

Despite the testing of the BabyCO that has been undertaken, normal operation of the BabyCO can be affected by other electrical/electronic equipment and portable and mobile RF communications equipment.

As the BabyCO is medical equipment, special precautions are needed regarding EMC (electromagnetic compatibility).

The BabyCo meter has an essential performance – the product should continue to operate as intended, to an accuracy of +/- 5% of reading or 1ppm (whichever is greater)"

It is important that the BabyCO is configured and installed/put into service, in accordance with the instructions/guidance provided herein and is used only in the configuration as supplied.

Changes or modifications to the BabyCO may result in increased emissions or decreased immunity of the BabyCO in relation to EMC performance.

The BabyCO should be used only with the PC serial cable provided by Carefusion (Cat. No. 36-CAB1000, Description: Interface cable for COBRA software). This cable should not be extended by the user. This cable should not be used with devices other than the BabyCO. If the cable is extended by the user, this may result in an increased level of emissions or decreased level of immunity, in relation to the BabyCO's EMC. Use of the cable with devices other than the BabyCO, may result in an increased level of emissions or decreased level of immunity, in relation to the other devices's EMC.

**WARNING:** The BabyCO should not be used adjacent to or stacked with other equipment. If adjacent or stacked use with other equipment is necessary, the BabyCO and the other equipment should be observed / monitored, to verify normal operation in the configuration in which it will be used.

<b>Guidance and Manufacturer's Declaration - Electromagnetic Emissions</b>			
<b>Emission Test</b>		<b>Compliance</b>	<b>Electromagnetic Environment - Guidance</b>
RF emissions CISPR 11		Group 1	The BabyCO uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11		Group B	The BabyCO is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC61000-3-2		Not Applicable (battery powered)	
Voltage fluctuations / flicker emissions IEC61000-3-3		Not Applicable (battery powered)	

<b>Guidance and Manufacturer's Declaration - Electromagnetic Immunity</b>				
<b>Immunity Test</b>		<b>IEC 60601 test level</b>	<b>Compliance Level</b>	<b>Electromagnetic Environment - Guidance</b>
Electrostatic discharge (ESD) IEC61000-4-2		± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC61000-4-4		± 2 kV for power supply lines  ± 1 kV for input / output lines	Not Applicable (Battery Powered and Serial cable <3m)	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC61000-4-5		± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not Applicable (battery Powered)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions		<5% U <sub>T</sub> (>95 % dip in U <sub>T</sub> ) For 0.5 cycle	Not Applicable (battery Powered)	Mains power quality should be that of a typical commercial or hospital

and voltage variations on power supply input lines IEC61000-4-11	40% $U_T$ (60 % dip in $U_T$ ) for 5 cycles 70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles <5% $U_T$ (>95 % dip in $U_T$ ) For 5 s		environment. If the user of the BabyCO requires continued operation during power mains interruptions, it is recommended that the BabyCO be powered from an uninterruptable power supply or a battery
Power frequency (50/60Hz) Magnetic field IEC61000-4-8	3 A/m	3 A/m	If incorrect operation occurs, it may be necessary to position the BabyCo further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.

NOTE  $U_T$  is the a.c. mains voltage prior to application of the test level.

#### Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The BabyCO is intended for use in the electromagnetic environment specified below. The customer or the user of the BabyCO should assure that it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the BabyCO, including any cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. <b>Recommended separation distance (<math>d</math>)</b> $d = 1.2\sqrt{P}$
Radiated RF IEC61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$

			<p>GHz</p> <p>Where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in metres (m).</p> <p>Fields strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<p><sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BabyCO is used exceeds the applicable RF compliance level above, the BabyCO should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orientating or relocating the BabyCO.</p> <p>b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

### **Recommended separation distances between portable and mobile RF communications equipment and the BabyCO**

The BabyCO is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the BabyCO can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BabyCO as recommended below, according to the maximum output power of the communications equipment.

<b>Rated Maximum Output Power of Transmitter in Watts (W)</b>	<b>Separation Distance in Meters (m) according to Frequency of Transmitter</b>
150 KHz to 80 MHz	80 MHz to 800 MHz      800 MHz to 2.5 GHz

	$d = 1.2\sqrt{P}$	$d = 1.2\sqrt{P}$	$d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures objects and people

## Symbols

Type B device

In accordance with Directive 93/42/EEC

0086

Disposal in compliance with WEEE

Consult the instructions for use

Caution: Consult the accompanying documents

Manufacturer

Date of Manufacture

**SN**

Serial number

Single patient use

Federal U.S. law restricts this device to sale by or on the order of a physician. (Rx only)

Medical Device listing mark for US and Canada by SGS Testing Services recognised by the American Occupational Safety and Health Administration (OSHA) for electrical safety and compliance.

## Specifications

Sensor type	Electro-chemical fuel cell
Range	0 - 100 ppm
Resolution	1 ppm
Green indicator light	0 to 6ppm
Amber indicator light	7 to 10ppm
Red indicator light	11 to 20ppm
Flashing red light	20+
Accuracy	+/-5% of reading or 1ppm whichever is the greater
Sensitivity drift	0.5%/ $^{\circ}$ C
Sensor life	2 to 5 years
Response time	< 15 sec (to 90% of reading)
Hydrogen cross sensitivity	<15%
Operating temperature	15 - 25 $^{\circ}$ C
Operating pressure	Atmospheric +/- 10%
Pressure coefficient	0.02% signal per mBar
Relative humidity	15 - 90% continuous
(Non condensing)	(0 - 99% intermittent)
Baseline drift	0ppm (auto-zero)
Long term drift	< 2% signal loss per month
Power source	Single Alkaline 9 volt PP3
Main battery life	30 hours of continuous use equivalent to approximately 2000 tests
Internal Battery	Lithium ½ AA 3.6 Volt
Internal battery life	10 years
Weight	180 g (Including battery)
Dimensions	170 x 60 x 26 mm
Display	3 ½ digit LCD
Storage & Transport temperature	-20 $^{\circ}$ to +70 $^{\circ}$ C
Storage & Transport humidity	30% to 90%

## **Consumables / Accessories**

<b>Cat. No.</b>	<b>Description</b>
36-PSA1200	22mm paediatric disposable mouthpieces (250 per box)
36-PSA2000	22mm mouthpiece adaptor with one-way valve (pack of 10)
36-PSA1100	Adapter for paediatric disposable mouthpieces
36-BAT1002	Alkaline PP3 Battery
36-MEC1184	Calibration Tool  Medican calibration gas canister
36-MCG020	20 litres of gas compressed in a 1-litre canister 20ppm Carbon Monoxide in air
36-MGA222	MicroFlow pressure reducer for Medican calibration gas canister with flow indicator
36-MEC1007	22mm reducing connector for calibration
36-CAB1000	Interface cable for COBRA software
36-SSC5000A	Perasafe sterilising powder (to make up 5 litres of solution)
36-VOL2104	Nose Clips (pack of 5)
V-861449	Silicone adapter "oval" for MicroGard® IIB
861427	Calibration syringe adapter

For further information or to place an order for Disposables/  
Supporting products please contact CAREFUSION, your local  
distributor or view our website [www.carefusion.com/micromedical](http://www.carefusion.com/micromedical)

**PLEASE NOTE: USE ONLY CAREFUSION ACCESSORIES**

**Please Note:** Information in this manual is subject to change without notice and does not represent commitment on the part of CareFusion 232 UK Ltd. The software may be used or copied only in accordance with the terms of that agreement. No part of this manual may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying and recording for any purpose without the written permission of CareFusion 232 UK Ltd.

## **Customer contact information**

### **UK Customers only**

For all Sales Order processing for products, training and spare parts, Service and Technical Support enquiries, please contact the following:

CareFusion UK 232 Ltd  
UK Customer Service & Support  
The Crescent  
Jays Close  
Basingstoke  
RG22 4BS

Customer Service Sales Enquiries:

Telephone: 01256 388550  
Email: [micro.uksales@carefusion.com](mailto:micro.uksales@carefusion.com)

Factory Repair and Administration Enquiries:

Telephone: 01256 388552  
Email: [micro.ukservice@carefusion.com](mailto:micro.ukservice@carefusion.com)

Technical Support Enquiries:

Telephone: 01256 388551  
Email: [support.rt.eu@carefusion.com](mailto:support.rt.eu@carefusion.com)

## **International customers only**

For all Sales Order processing for products and Spare parts, Service and Technical Support enquiries, please contact the following:

CareFusion Germany 234 GmbH  
Customer Service & Support International  
Leibnizstraße 7  
97204 Hoechberg  
Germany

Customer Service Sales Enquiries:

Telephone: 0049 931 4972 670  
Email: [micro.internationalsales@carefusion.com](mailto:micro.internationalsales@carefusion.com)

Factory Repair and Administration Enquiries:

Telephone: 0049 931 4972 867  
Email: [support.admin.eu@carefusion.com](mailto:support.admin.eu@carefusion.com)

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## **Quellen**

3. **Carbon monoxide in breath in relation to smoking and carboxyhemoglobin levels.**  
Wald NJ, Idle M, Boreham J, Baily A  
Thorax 1981; 36, 366-369
4. **Definition of a reliable threshold value for detecting current smokers by CO measurement**  
Marino Luigi; Latini Roberto; Barbano Gina; Bazzlerla Giorgio;  
De Luca Anita, Nardini Stefano – Respiratory and TB Unit-  
General Hospital- Via forlanini, 71-I-31029-Vittorio Veneto (TV-  
ITALY).

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## Referencias

3. Gomez C, Berlin I, Marquis P, Delcroix M.  
**Expired air carbon monoxide concentration in mothers and their spouses above 5 ppm is associated with decreased foetal growth.**  
Preventative Medicine 40 (2005) 10-15

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## Referencias

5. **Carbon monoxide in breath in relation to smoking and carboxyhemoglobin levels.**  
Wald NJ, Idle M, Boreham J, Baily A  
Thorax 1981; 36, 366-369
6. **Definition of a reliable threshold value for detecting current smokers by CO measurement**  
Marino Luigi; Latini Roberto; Barbano Gina; Bazzerla Giorgio;  
De Luca Anita, Nardini Stefano – Respiratory and TB Unit-  
General Hospital- Via forlanini, 71-I-31029-Vittorio Veneto (TV-  
ITALY).

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## Referenties

9. **Carbon monoxide in breath in relation to smoking and carboxyhemoglobin levels.**  
Wald NJ, Idle M, Boreham J, Baily A  
Thorax 1981; 36, 366-369
10. **Definition of a reliable threshold value for detecting current smokers by CO measurement**  
Marino Luigi; Latini Roberto; Barbano Gina; Bazzlerla Giorgio;  
De Luca Anita, Nardini Stefano – Respiratory and TB Unit-  
General Hospital- Via forlanini, 71-I-31029-Vittorio Veneto (TV-  
ITALY).

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## Riferimenti

11. **Carbon monoxide in breath in relation to smoking and carboxyhemoglobin levels.**  
Wald NJ, Idle M, Boreham J, Baily A  
*Thorax* 1981; 36, 366-369
12. **Definition of a reliable threshold value for detecting current smokers by CO measurement**  
Marino Luigi; Latini Roberto; Barbano Gina; Bazzlerla Giorgio;  
De Luca Anita, Nardini Stefano – Respiratory and TB Unit-  
General Hospital- Via forlanini, 71-I-31029-Vittorio Veneto (TV-  
ITALY).

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© CareFusion 2013