CareFusion

MicroRPM

(Respiratory Pressure Meter)

Operating Manual – English

Operating Manual – U.S

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Contents

I	Introduction - English	4
	Contents	5
	PUMA PC Software	6
	Contraindications	7
	Warnings and Cautions	7
	Indications for Use	8
	Operation – Mouth Pressures (PI _{max} /MIP + PE _{max} /MEP)	9
	PI _{max} (MIP) Test	9
	PE _{max} (MEP) Test	10
	Operation – SNIP (Sniff Nasal Inspiratory Pressure)	11
	SNIP Test	12
	Switching Off	12
	Battery	13
	Battery Replacement	13
	Cleaning Instructions	14
	External Surfaces of the Spirometer	14
	Cleaning accessories	14
	Calibration	15
	Servicing	16
	Trouble Shooting Information	17
	Safety Designation per IEC 60601-1	18
	Electromagnetic Compatibility (EMC) to EN60601-1:2007	19
	Symbols	24
	Specifications	25
	Consumables and Accessories	26
II.	Introduction – U.S.	29
	Package Contents	30
	PUMA PC Software	31
	Contraindications	32
	Warnings and Cautions	32

Indications for Use	33
Operation – Mouth Pressures.....	34
Operation – Sniff Nasal Inspiratory Pressure	36
Switching Off	37
Battery Low Voltage Indication	37
Battery Replacement.....	38
Cleaning Instructions	39
External Surfaces of the Spirometer	39
Cleaning accessories	39
Calibration	40
Servicing	42
Trouble Shooting Information.....	42
Safety Designation per IEC 60601-1	43
Electromagnetic Compatibility (EMC) to EN60601-1:2007	44
Symbols	49
Specifications	50
Consumables / Supporting Products	51
III. Introduction - Français	52
Contenu De L’emballage	53
Le logiciel PUMA	55
Contre-indications.....	56
Avertissements et mises en garde.....	56
Indications d'utilisation.....	57
Mode D’emploi – Pressions Buccales.....	58
Mode D’emploi – Pression Inspiratoire Sniff Nasale	60
Mise hors tension.....	62
Témoin Du Niveau De Charge De La Pile.....	62
Remplacement de la pile.....	62
Instructions de nettoyage	63
Surfaces extérieures du spiromètre.....	64
Accessoires de nettoyage.....	64
Etalonnage	65

Entretien	67
Diagnostic des pannes	67
Certificat de sécurité selon la norme CEI 60601-1	68
Compatibilité électromagnétique (CEM) à la norme EN60601-1:2007	69
Symboles.....	74
Spécifications	76
Pièces De Rechange / Accessoires.....	77
IV Einführung - Deutsch.....	79
Inhalt.....	80
PUMA-PC-Software.....	81
Gegenanzeigen.....	82
Warn- und Vorsichtshinweise	82
Gebrauchsanweisung	84
Betrieb – Munddruck (P _I max/MIP + P _E max/MEP)	85
P _I max (MIP) Test.....	86
P _E max (MEP) Test	87
Betrieb – SNIP (Einatemnasendruck).....	88
SNIP Test	89
Ausschalten	89
Batterie	90
Batteriewechsel	90
Reinigungsanleitung.....	91
Außenoberflächen des Spirometers	91
Reinigungszubehör.....	91
Kalibrierung	92
Wartung	94
Informationen zur Fehlerbehebung	94
Sicherheitsbezeichnung gemäß IEC 60601-1	95
Elektromagnetische Verträglichkeit (EMV) gemäß EN60601-1:2007	96
Symbole.....	103

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Operating manual

I Introduction - English

The respiratory pressure meter is a hand held instrument designed for rapid assessment of inspiratory and expiratory muscle strength. The unit can measure the maximum inspiratory and expiratory mouth pressures, MIP and MEP, and the Sniff Nasal Inspiratory Pressure, SNIP. The result of each measurement is presented in units of cmH₂O gauge pressure on the liquid crystal display screen.

The unit is easy to operate, battery powered, and is supplied with all the necessary attachments required for immediate use.


The functionality of the unit may be greatly increased when connected to a PC running PUMA software. This application has many advanced features including:

- Real time display of pressure/time curves
- Overlay of successive curves
- Predicted values
- Patient database
- Incentive display
- Manoeuvre quality check
- Manoeuvre variability measurement


Contents

The MicroRPM is supplied as follows:

- 1. MicroRPM Microcomputer Unit
- 2. Rubber Flanged Mouthpiece Cat. No. 36-MTH6400
- 3. Alkaline battery (type 6LR61) Cat. No. 36-BAT1002
- 4. Expiratory Pressure Valve Assembly Cat. No. 36-ASS1221
- 5. Inspiratory Pressure Valve Assembly Cat. No. 36-ASS1222
- 6. Calibration Screwdriver Cat. No. 36-MEC1299
- 7. Nasal Probes Size 1-4
- 8. Nasal Probe Adapter Cat. No. 36-ASS1091
- 9. Mouth Pressure Bacterial Filters Cat. No. 36-FIL6050 (50)

PUMA PC Software

The functionality of the portable MicroRPM is greatly increased when connected to the PUMA PC Software, via an RS232 cable to the Serial Port on the side of the unit.

PUMA PC Software is available as a free download from www.carefusion.com/micromedical. Please contact Technical Services for installation instructions and system requirements.

The PUMA PC Software is a unique, user friendly multi windows platform for the performance, storage and analysis of the respiratory muscle strength measurements of PImax or MIP (Maximum Inspiratory Pressure), PEmax or MEP (Maximum Expiratory Pressure) and SNIP (Sniff Nasal Inspiratory Pressure). In addition, PUMA PC Software offers the user features such as live graphical displays, predicted values, printing formats, incentives, trends, post medication or exercise comparisons and additional fatigue indicators.

Note: The Respiratory Pressure Meter should only be connected to a computer that is manufactured in accordance with EN 60601-1.

Note: Keep the PC out of reach of the patient at all times.

Contraindications

- Pathological conditions resulting in relatively large pressure swings in the abdomen or thorax
- Aneurisms
- Uncontrolled hypertension
- Urinary incontinence

Warnings and Cautions

The following terms are used as follows in this manual:

Caution: Possibility of injury or serious damage

Warning: conditions or practices that could result in personal injury

Please Note: Important information for avoiding damage to the instrument or facilitating operation of the instrument.

Note: The device should only be used by qualified personnel trained in lung function testing.

CAUTION: Read the manual before use

CAUTION: For batteries, do not attempt to charge, connect improperly, or dispose of in fire as there is possibility of leakage or explosion. Follow manufacturer's recommendation for proper disposal.

WARNING: The instrument is not suitable for use in the presence of explosive or flammable gases, flammable anaesthetic mixtures or in oxygen rich environments.

CAUTION: Bacterial filters are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may increase air resistance and lead to an incorrect measurement.

PLEASE NOTE: The product and the battery you have purchased should not be disposed of as unsorted waste. Please utilise your local WEEE collection facilities for the disposal of this product.

PLEASE NOTE: Degree of protection against Ingress of Water is IPX0.

CAUTION: When you connect the Respiratory Pressure Meter to other equipment, always make sure that the whole combination complies with the international safety standard IEC 60601-1 for medical electrical systems. During measurements, connect the Respiratory Pressure Meter only to computers that comply with IEC/EN 60601-1 / ANSI/AAMI ES60601-1:2005.


Indications for Use

The Micro RPM (Respiratory Pressure Meter) is a hand held diagnostic instrument designed for rapid assessment of inspiratory and expiratory muscle strength. The unit can measure the maximum inspiratory and expiratory mouth pressures, MIP and MEP, and the Sniff Nasal inspiratory Pressure, SNIP. The system is intended for use with pediatric and adult patients over the age of 3 years in hospitals, physician offices, laboratories and occupational health testing environments.

Operation – Mouth Pressures (P_Imax/MIP + P_Emax/MEP)

Insert the battery into the compartment at the rear of the MicroRPM.

Fit the required Pressure Valve Assembly ('Inspiratory' for P_Imax (MIP), 'Expiratory' for P_Emax (MEP)) into the MicroRPM; insert a new Bacterial Filter into the Pressure Valve Assembly and then the Rubber Flanged Mouthpiece onto the Bacterial Filter, as shown below.

P_Imax (MIP) Test

Slide the MicroRPM switch from 'Off' to 'MIP/MEP', whilst applying no pressure to the mouthpiece. Rotating segments will be displayed whilst the unit performs an auto-zero function.

When the MicroRPM is ready a 'beep' will be heard and '0' displayed.

To perform the test instruct the subject to insert the mouthpiece into the mouth, ensuring the flange is positioned over the gums and inside the lips, whilst the 'bite blocks' are between the teeth.

The subject should then exhale to RV (Residual Volume), lungs empty, then perform a 'Mueller' manoeuvre, a forced inhalation against the MicroRPM with as much effort as possible for as long as possible (minimum 2 seconds).

The display will report the result, the maximum average inspiratory pressure sustained over a 1 second period of the test, in centimetres of water (cmH₂O). Ideally, the subject should repeat this test 3 times to ascertain a best value.

PE_{max} (MEP) Test

Slide the MicroRPM switch from 'Off' to 'MIP/MEP', whilst applying no pressure to the mouthpiece. Rotating segments will be displayed whilst the unit performs an auto-zero function.

When the MicroRPM is ready a 'beep' will be heard and '0' displayed.

To perform the test instruct the subject to insert the mouthpiece into the mouth, ensuring the flange is positioned over the gums and inside the lips, whilst the 'bite blocks' are between the teeth.

The subject should then inhale to TLC (Total Lung Capacity), lungs full, then perform a 'Valsalva' manoeuvre, a forced exhalation against the MicroRPM with as much effort as possible for as long as possible (minimum 2 seconds).

The display will report the result, the maximum average expiratory pressure sustained over a 1 second period of the test, in cmH₂O. Ideally, the subject should repeat this test 3 times to ascertain a best value.


To repeat either the PImax or PEmax tests the MicroRPM must firstly be returned to the 'Off' position.

Operation – SNIP (Sniff Nasal Inspiratory Pressure)

Insert the battery into the rear of the MicroRPM.

Fit the Nasal Probe Adapter into the MicroRPM and then attach the correct size Nasal Probe, as shown below.

The correct size (1-4) can be ascertained by fitting a Nasal Probe to the unit, then firmly inserting the Nasal Probe into a nostril. Instruct the subject to block the open nostril with a finger and then to attempt a sniff. The correct Nasal Probe size has been selected once there is no leakage around the Nasal Probe.

SNIP Test

Slide the MicroRPM switch from 'Off' to 'SNIP', whilst applying no pressure to the Nasal Probe. Rotating segments will be displayed whilst the unit performs an auto-zero function.

When the MicroRPM is ready a 'beep' will be heard and '0' displayed.

To perform the test, instruct the subject to insert the chosen Nasal Probe firmly into a nostril, whilst ensuring the other nostril remains open throughout the test.

The subject should then breathe normally and at the end of a normal tidal expiration, FRC (Functional Residual Capacity), be instructed to perform a short, sharp voluntary sniffing manoeuvre with as much effort as possible.

The display will report the result, the peak inspiratory nasal pressure, in cmH₂O.

On subsequent tests the MicroRPM will continue to display the highest SNIP value, overwriting previous values. Ideally, the subject should repeat this test 10-15 times to ascertain the highest value.

Switching Off

The MicroRPM is switched off by sliding the switch back to the 'Off' position.

Battery

The battery level is automatically checked when the unit is switched on.

When the battery is low 'bAt' will be displayed before the auto-zero procedure. The MicroRPM may be used when this occurs, but the battery should be replaced as soon as possible.

When the battery is completely exhausted the unit will beep twice and turn itself off immediately.

Battery Replacement

Locate the sliding cover situated on the rear of the unit, towards the bottom of the device.

Place your thumb over the round thumb indent, press gently and slide the cover to the right to remove it from the unit.

Lift the old battery out, and holding the battery terminal by the plastic body, pull it off the old battery.

Plug the new battery into the battery terminal, taking care that the correct polarity is observed.

Push the battery back into the battery holder, and replace the battery cover onto the guides. Slide the battery cover to the left until it is fully home.

Note: Please remove the battery if the meter is likely to be unused for some time.

CAUTION: Do not open the battery cover when the device is turned on.

CAUTION: The operator should not touch the contacts of the battery and the patient at the same time.

Please Note: Dispose of the waste battery in accordance with EU Waste Battery Regulations.

Cleaning Instructions

Disinfection of contaminated parts is only effective after having them carefully pre-cleaned. CareFusion recommends the tested solution of PeraSafe sterilizing powder (36-SSC5000A) for pre-cleaning and disinfection. If a different solution is used, please follow the given manufacturer's instructions.

The device must not be wiped with any aqueous solutions and must not be exposed to solvents e.g. alcohol, chloride solutions as there are electronic components inside that will be permanently damaged.

CAUTION: Switch off the device before cleaning.

External Surfaces of the Spirometer

CAUTION: Do not attempt to wash or immerse the Respiratory Pressure Meter in water or cleaning fluid, as there are electronic components inside that will be permanently damaged

The external housing of the pressure meter may be wiped with sterile wipes or a damp cloth that has been immersed in a cold sterilising solution.

Cleaning accessories

The MicroRPM unit is protected from contamination by the Bacterial Filter (36-FIL6050) during mouth pressures measurements.

The Rubber Flanged Mouthpiece (36-MTH6400), the Expiratory and Inspiratory Pressure Valve Assembly (36-ASS1221, 36-ASS1222) and the Nasal Probes (36-NPROBE01, 36-NPROBE02, 36-NPROBE03, 36-NPROBE04), however, may be immersed in a cold sterilising solution such as PeraSafe (36-SSC5000A). Rinse thoroughly and leave to dry before reassembly.


Important note: Used mouthpieces and Nasal Probes, which are not sterilized, must be disposed off immediately after each use.

If there are changes on the material surfaces (cracks, brittleness) the respective parts must be disposed off.

Calibration

The calibration is factory set and should remain stable indefinitely.


The calibration may be checked and altered, however, by connecting the device via the Pressure Calibration Kit (36-ASS1234) to a Water Manometer, as shown below:

Slide the MicroRPM switch from 'Off' to 'SNIP' whilst applying no pressure to the unit.

Gently fill the syringe until negative pressure of between 200 and 300cmH₂O is read on the Water Manometer. Check that the reading on the MicroRPM is within 3% of this reading.

Please Note: Calibration can only be adjusted in a positive direction, as the MicroRPM measures the highest peak value. Therefore, if the reading was greater than that on the Water Manometer, the calibration screw must be turned anti-clockwise first, before calibration is attempted.

To calibrate the MicroRPM, connect to the Water Manometer as before, filling the syringe to obtain the required negative pressure. To change the calibration, slowly turn the calibration screw in a clockwise direction until the same value is displayed on the MicroRPM.

Servicing

If your unit requires service or repair please see page 137 for contact details.

A full service manual including circuit diagrams and parts list is available on request.

Trouble Shooting Information

Should you encounter problems operating your MicroRPM unit, please consult the table below:

Problem	Possible Cause	Solution
unit will not turn on	Battery is exhausted	Replace the battery
	Slide switch connection	Return unit for servicing.
Display shows reading before test has been performed	Internal tubing to pressure sensor kinked Return unit for servicing	

Electromagnetic Compatibility (EMC) to EN60601-1:2007

WARNING: Use of portable phones or other radio frequency (RF) emitting equipment near the system may cause unexpected or adverse operation

The MicroRPM has been tested to EN60601-1-2:2007, regarding its ability to operate in an environment containing other electrical/electronic equipment (including other medical devices).

The purpose of this testing is to ensure that the MicroRPM is not likely to adversely affect the normal operation of other such equipment and that other such equipment is not likely to adversely affect the normal operation of the MicroRPM.

Despite the testing of the MicroRPM that has been undertaken, normal operation of the MicroRPM can be affected by other electrical/electronic equipment and portable and mobile RF communications equipment.

As the MicroRPM is medical equipment, special precautions are needed regarding EMC (electromagnetic compatibility).

It is important that the MicroRPM is configured and installed/put into service, in accordance with the instructions/guidance provided herein and is used only in the configuration as supplied.

Changes or modifications to the MicroRPM may result in increased emissions or decreased immunity of the MicroRPM in relation to EMC performance.

The MicroRPM should be used only with the accessories (RS232 cables) supplied (which are referenced in the accessories section of this manual). None of the MicroRPM cables should be extended in length by the user.

If any cables are extended by the user or non approved accessories are used, this may result in an increased level of emissions or decreased level of immunity, in relation to the MicroRPM EMC. None of the MicroRPM accessories should be used with other devices, as this may result in an increased level of emissions or decreased level of immunity, in relation to the other devices' EMC.

The MicroRPM has an essential performance – The respiratory pressure readings on the product must remain within a tolerance of +/- 3% and the unit firmware must not cease responding. Warning: In the event the product is operated in the presence of significant Electromagnetic Fields (particularly in the frequency range 40-60MHz), whilst in the PC connected mode, ensure the results on the unit and PC are the same. If the results differ then re-locating the product away from sources of interference should resolve any issue.

WARNING: The MicroRPM should not be used adjacent to or stacked with other equipment. If adjacent or stacked use with other equipment is necessary, the MicroRPM and the other equipment should be observed / monitored, to verify normal operation in the configuration in which it will be used.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions			
The MicroRPM is intended for use in the electromagnetic environment specified below. The customer or the user of the MicroRPM should assure that it is used in such an environment			
Emission Test	Compliance	Electromagnetic Environment - Guidance	
RF emissions CISPR11	Group 1	The MicroRPM uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF emissions CISPR11	Class B	The MicroRPM is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purpose.	
Harmonic Emissions IEC 61000-3-2	[Not Applicable]		
Voltage fluctuations / flicker emissions IEC 61000-3-3	[Not applicable]		

Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
The MicroRPM is intended for use in the electromagnetic environment specified below. The customer or the user of the MicroRPM should assure that it is used in such an environment.			
Immunity Test	IEC 60601 test level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6 kV contact	+/- 6 kV contact	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
	+/- 8 kV air	+/- 8 kV air	
Electrical fast transient/burst IEC 61000-4-4	+/- 2 kV for power supply lines	+/- 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
	+/- 1 kV for input/output lines		

Surge IEC 61000-4-5	+/- 1 kV line(s) to line (s) +/- 2 kV line(s) to earth	Not Applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% UT (>95% dip in UT) for 0.5% cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (>95% dip in UT) for 5s	Not Applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the micro RPM requires continued operation during power main interruptions, it is recommended that the micro RPM be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristics of a typical location in a typical commercial or hospital environment.
NOTE UT is the a.c. mains voltage prior to application of the test level.			

Guidance and Manufacturer’s Declaration – Electromagnetic Immunity

The MicroRPM is intended for use in the electromagnetic environment specified below. The customer or the user of the MicroRPM should assure that it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance Level	Electromagnetic Environment -Guidance
<p>Conducted RF IEC61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150kHz to 80MHz</p> <p>3 V/m 80MHz to 2.5GHz</p>	<p>3 V</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the micro RPM, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p>$d = 1.2 \times \sqrt{P}$</p> <p>$d = 1.2 \times \sqrt{P}$ 80MHz to 800MHz</p> <p>$d = 2.3 \times \sqrt{P}$ 800MHz to 2.5GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from a fixed RF transmitters as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.
 NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radio, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the MicroRPM is used exceeds the applicable RF compliance level above, the MicroRPM should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orientating or relocating the MicroRPM.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m

Recommended separation distances between portable and mobile RF communications equipment and the MicroRPM

The MicroRPM is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the MicroRPM can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the MicroRPM as recommended below, according to the maximum output power of the communications equipment

Rated Maximum Output Power of Transmitter in Watts (W)	Separation Distance in Meters (m) according to Frequency of Transmitter		
	150 KHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1.2 \times \sqrt{P}$	$d = 1.2 \times \sqrt{P}$	$d = 2.4 \times \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures objects and people.

Symbols

Type B device

0086

In accordance with Directive 93/42/EEC

Disposal in compliance with WEEE

Consult the instructions for use

Caution: Consult the accompanying documents

Date of Manufacture

Manufacturer

Serial number

Single patient use

Federal U.S. law restricts this device to sale by or on the order of a physician. (Rx only)

Nationally Recognized Test Laboratory (NRTL) officially recognized by the American Occupational Safety and Health Administration (OSHA) for OSHA Electrical safety and compliance.

Specifications

Measurements	Maximum Expiratory Pressure (MEP) Maximum Inspiratory Pressure (MIP) Sniff Nasal Inspiratory Pressure (SNIP)
With Puma™	Maximum Rate of Pressure Development (MRPD) Maximum Rate of Relaxation (MRR) Tau
Operating Pressure:	+/- 300cmH2O (±5PSID)
Burst Pressure:	+/- 700cmH2O (±20PSID)
Accuracy:	+/-3%
Resolution:	1cmH2O
Power Supply:	Single Alkaline 9V 6LR61
Dimensions:	170 x 60 x 26mm
Weight (unit with battery):	175g
Weight (in carry case) :	750g
Operating Temperature :	0 - 40 °C
Operating Humidity:	30% - 90% RH
Storage & Transport temperature:	-20 - +70 ° C
Storage & Transport humidity:	10% - 90% RH

Consumables and Accessories

Cat. No.	Description
36-FIL6050	Bacterial Filters (50 per box)
36-MTH6400	Rubber Flanged Mouthpiece
36-ASS1221	Expiratory Pressure Valve Assembly
36-ASS1222	Inspiratory Pressure Valve Assembly
36-NPROBE01	Extra Small Nasal Probes (10 per box)
36-NPROBE02	Small Nasal Probes (10 per box)
36-NPROBE03	Medium Nasal Probes (10 per box)
36-NPROBE04	Large Nasal Probes (10 per box)
36-ASS1091	Nasal Probe Adapter
36-MEC1299	Calibration Screwdriver
36-BAT1002	Alkaline Battery (type 6LR61)
36-SSC5000A	PeraSafe Sterilising Powder – 81g (Makes 5 litres of solution)
36-ASS1234	Pressure Calibration Kit
36-CAB1000	Interface Cable for PUMA Software

For further information or to place an order for Disposables/
Supporting products please contact CAREFUSION, your local
distributor or view our website www.carefusion.com/micromedical

PLEASE NOTE: USE ONLY CAREFUSION ACCESSORIES

Please Note: Information in this manual is subject to change without notice and does not represent commitment on the part of CareFusion 232 UK Ltd. The software may be used or copied only in accordance with the terms of that agreement. No part of this manual may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying and recording for any purpose without the written permission of CareFusion 232 UK Ltd.

MicroRPM

(Respiratory Pressure Meter)

Operating manual

Federal (USA) law restricts this device to sale by or on the order of a physician or licensed practitioner.

CareFusion,
22745 Savi Ranch Parkway,
Yorba Linda
CA 92887-4668
USA

II. Introduction – U.S

The respiratory pressure meter is a hand held instrument designed for rapid assessment of inspiratory and expiratory muscle strength. The unit can measure the maximum inspiratory and expiratory mouth pressures, MIP and MEP, and the Sniff Nasal Inspiratory Pressure, SNIP. The result of each measurement is presented in units of cmH₂O gauge pressure on the liquid crystal display screen.

The unit is easy to operate, battery powered, and is supplied with all the necessary attachments required for immediate use.


The functionality of the unit may be greatly increased when connected to a PC running PUMA software. This application has many advanced features including:

- Real time display of pressure/time curves
- Overlay of successive curves
- Predicted values
- Patient database
- Incentive display
- Manoeuvre quality check
- Manoeuvre variability measurement


Package Contents

The respiratory pressure meter is supplied with the following items:-

1. Microcomputer unit
2. Rubber mouthpiece
3. Alkaline PP3 battery
4. Expiration pressure valve assembly
5. Inspiration pressure valve assembly
6. Calibration screwdriver
7. Nasal probe sizes 1 to 4
8. Nasal probe adapter
9. Mouth pressure bacterial filters

PUMA PC Software

The functionality of the portable MicroRPM is greatly increased when connected to the PUMA PC Software, via an RS232 cable to the Serial Port on the side of the unit.

PUMA PC Software is available as a free download from www.carefusion.com/micromedical. Please contact Technical Services for installation instructions and system requirements.

The PUMA PC Software is a unique, user friendly multi windows platform for the performance, storage and analysis of the respiratory muscle strength measurements of PImax or MIP (Maximum Inspiratory Pressure), PEmax or MEP (Maximum Expiratory Pressure) and SNIP (Sniff Nasal Inspiratory Pressure). In addition, PUMA PC Software offers the user features such as live graphical displays, predicted values, printing formats, incentives, trends, post medication or exercise comparisons and additional fatigue indicators.

Note: The Respiratory Pressure Meter should only be connected to a computer that is manufactured in accordance with EN 60601-1.

Note: Keep the PC out of reach of the patient at all times.

Contraindications

- Pathological conditions resulting in relatively large pressure swings in the abdomen or thorax
- Aneurisms
- Uncontrolled hypertension
- Urinary incontinence

Warnings and Cautions

The following terms are used as follows in this manual:

Caution: Possibility of injury or serious damage

Warning: conditions or practices that could result in personal injury

Please Note: Important information for avoiding damage to the instrument or facilitating operation of the instrument.

Note: The device should only be used by qualified personnel trained in lung function testing.

CAUTION: Read the manual before use

CAUTION: For batteries, do not attempt to charge, connect improperly, or dispose of in fire as there is possibility of leakage or explosion. Follow manufacturer's recommendation for proper disposal.

WARNING: The instrument is not suitable for use in the presence of explosive or flammable gases, flammable anaesthetic mixtures or in oxygen rich environments.

CAUTION: Bacterial filters are single patient use. If used on more than one patient, there is a risk of cross-infection. Repeated use may increase air resistance and lead to an incorrect measurement.

PLEASE NOTE: The product and the battery you have purchased should not be disposed of as unsorted waste. Please utilize your local WEEE collection facilities for the disposal of this product.

PLEASE NOTE: Degree of protection against penetration of water is IPX0.


CAUTION: When you connect the Respiratory Pressure Meter to other equipment, always make sure that the whole combination complies with the international safety standard IEC 60601-1 for medical electrical systems. During measurements, connect the Respiratory Pressure Meter only to computers that comply with IEC/EN 60601-1 / ANSI/AAMI ES60601-1:2005.

Indications for Use

The Micro RPM (Respiratory Pressure Meter) is a hand held diagnostic instrument designed for rapid assessment of inspiratory and expiratory muscle strength. The unit can measure the maximum inspiratory and expiratory mouth pressures, MIP and MEP, and the Sniff Nasal inspiratory Pressure, SNIP. The system is intended for use with adults and pediatric patients over the age of 3 years in hospitals, physician offices, laboratories and occupational health testing environments.

Operation – Mouth Pressures

Insert the battery into the battery compartment at the rear of the unit, fit the flanged rubber mouthpiece over the required valve assembly and insert into the respiratory pressure meter as shown below:

Two valve assemblies are provided. The one marked expiration pressure allows the patient to inspire fully through the valve, which then closes during expiration to allow the meter to measure the maximum expired pressure averaged over one second.

The valve marked inspiratory pressure works in the exact opposite manner.

To measure the maximum expiratory pressure (PE_{max}), slide the switch from the "off" position to the MIP/MEP position. Rotating segments will be displayed whilst the unit performs an auto-zero function. During this time there should be no applied pressure.

When the unit is ready for use 0 will be displayed.

To perform an expiratory pressure test instruct the subject to insert the mouthpiece into the mouth ensuring that the flange is positioned over the gums and inside the lips and that the "bite blocks" are between the teeth. They should then inhale to TLC (total lung capacity) and then exhale with as much effort as possible through the controlled leak of the meter for at least 2 seconds. The reading displayed is the maximum average expiratory pressure over 1 second.

For an accurate measurement it is important that there are no leaks in the system other than the small controlled leak of the instrument. This leak prevents the generation of a false high reading generated by closure of the glottis and compression of air in the mouth using the facial muscles.


Before repeating the measurement, or before performing an inspiratory manoeuvre, the instrument must be reset by sliding the switch to the "off" position.

To measure the maximum inspiratory pressure (PI_{max}) follow the above procedure but with the valve marked inspiratory pressure connected. Slide the switch on the face of the meter from 'off' position to the MIP/MEP position and instruct the subject to exhale to RV (residual volume) through the mouthpiece and then to inhale with as much effort as possible for at least 2 seconds.

The reading displayed is the maximum average inspiratory pressure that was sustained over 1 second.

Operation – Sniff Nasal Inspiratory Pressure

Before commencing a test the correct size of nasal probe to be used must be ascertained. The unit is supplied with four sizes designed to cover a wide range of nostrils. Select the size that appears to be the most appropriate, connect it to the unit as shown below, and ask the subject to insert the probe firmly into a nostril. Once fitted the subject should close the other nostril with a finger and try to inhale through the probe to check the sealing. When the correct size has been selected proceed as follows: Insert the battery into the battery compartment at the rear of the unit, fit the selected probe as shown below:

To measure the maximum sniff inspiratory pressure, slide the switch from the "off" position to the SNIP position. Rotating segments will be displayed whilst the unit performs an auto-zero function.

During this time there should be no applied pressure.

When the unit is ready for use 0 will be displayed.

To perform a SNIP test instruct the subject to insert the nasal probe into a nostril such that a good seal is made. After a normal expiration they should inhale with as much effort as possible through the free nostril with the mouth closed. The manoeuvre is performed from functional residual volume (FRC). The reading displayed is the peak inspiratory nasal pressure.

Before repeating the measurement, the instrument must be reset by sliding the switch to the "off" position.

Switching Off

The MicroRPM is switched off by sliding the switch back to the 'Off' position.

Battery Low Voltage Indication

The battery level is checked when the unit is switched on.

When the battery is nearly exhausted 'bAt' is displayed before the auto-zero procedure begins. The unit may be used when this occurs provided the test is performed immediately. The battery should be replaced as soon as possible.

When the battery is completely exhausted the unit will beep twice and turn itself off immediately after turning on.

Note: If the respiratory pressure meter is not used for long periods, the battery should be removed to prevent damage to the instrument by possible leakage.

Battery Replacement

Locate the sliding cover situated on the rear of the unit, towards the bottom of the device.

Place your thumb over the round thumb indent, press gently and slide the cover to the right to remove it from the unit.

Lift the old battery out, and holding the battery terminal by the plastic body, pull it off the old battery.

Plug the new battery into the battery terminal, taking care that the correct polarity is observed.

Push the battery back into the battery holder, and replace the battery cover onto the guides. Slide the battery cover to the left until it is fully in place.

Note: Please remove the battery if the meter is likely to be unused for some time.

CAUTION: Do not open the battery cover when the device is turned on.

CAUTION: The operator should not touch the contacts of the battery and the patient at the same time.

Please Note: Dispose of the waste battery in accordance with EU Waste Battery Regulations.

Cleaning Instructions

Disinfection of contaminated parts is only effective after careful preliminary cleaning. CareFusion recommends the tested solution of PeraSafe sterilizing powder (36-SSC5000A) for pre-cleaning and disinfection. If a different solution is used, please follow the given manufacturer's instructions.

The device must not be wiped with any aqueous solutions and must not be exposed to solvents e.g. alcohol, chloride solutions, as there are electronic components inside that will be permanently damaged.

CAUTION: Switch off the device before cleaning.

External Surfaces of the Spirometer

CAUTION: Do not attempt to wash or immerse the Respiratory Pressure Meter in water or cleaning fluid, as there are electronic components inside that will be permanently damaged

The external housing of the pressure meter may be wiped with sterile wipes or a damp cloth that has been immersed in a cold sterilizing solution.

Cleaning accessories

The MicroRPM unit is protected from contamination by the Bacterial Filter (36-FIL6050) during mouth pressures measurements.

The Rubber Flanged Mouthpiece (36-MTH6400), the Expiratory and Inspiratory Pressure Valve Assembly (36-ASS1221, 36-ASS1222) and the Nasal Probes (36-NPROBE01, 36-NPROBE02, 36-NPROBE03, 36-NPROBE04), however, may be immersed in a cold sterilizing solution such as PeraSafe (36-SSC5000A). Rinse thoroughly and allow to dry before reassembly.


Important note: Used mouthpieces and Nasal Probes, which are not sterilized, must be disposed off immediately after each use.

If there are changes on the material surfaces (cracks, brittleness) the respective parts must be disposed off.

Calibration

Calibration is factory set and should remain stable indefinitely.

However, calibration may be checked by connecting the device to a manometer as shown below:


Turn the respiratory pressure meter on to the SNIP position.

Very slowly fill the syringe until a negative pressure of between 200 and 300cm is obtained on the manometer.

Check that the reading on the respiratory pressure meter is within 3% of the manometer reading.

If adjustment is required then the following procedure must be followed.

The calibration may only be adjusted in a positive direction as the meter monitors the peak pressure value. Therefore if the respiratory pressure meter reading was greater than the manometer reading the calibration screw must be turned anti-clockwise a few turns before calibration is attempted.

Connect the respiratory pressure meter to the manometer as previously shown. Fill the syringe to obtain the required negative pressure and then turn the calibration screw slowly in a clockwise direction until the same value is displayed on the meter.

Servicing

If your unit requires service or repair please see page 137 for contact details.

A full service manual including circuit diagrams and parts list is available on request.

Trouble Shooting Information

Should you encounter problems operating your MicroRPM unit, please consult the table below:

Problem	Possible Cause	Solution
unit will not turn on	Battery is spent/discharged	Replace the battery
	Slide switch connection	Return unit for servicing.
Display shows reading before test has been performed	Internal tubing to pressure sensor kinked Return unit for servicing	

Electromagnetic Compatibility (EMC) to EN60601-1:2007

WARNING: Use of portable phones or other radio frequency (RF) emitting equipment near the system may cause unexpected or adverse operation

The MicroRPM has been tested to EN60601-1-2:2007, regarding its ability to operate in an environment containing other electrical/electronic equipment (including other medical devices).

The purpose of this testing is to ensure that the MicroRPM is not likely to adversely affect the normal operation of other such equipment and that other such equipment is not likely to adversely affect the normal operation of the MicroRPM.

Despite the testing of the MicroRPM that has been performed, normal operation of the MicroRPM can be affected by other electrical/electronic equipment and portable and mobile RF communications equipment.

As the MicroRPM is medical equipment, special precautions are needed regarding EMC (electromagnetic compatibility).

It is important that the MicroRPM is configured and installed/put into service, in accordance with the instructions/guidance provided herein and is used only in the configuration as supplied.

Changes or modifications to the MicroRPM may result in increased emissions or decreased immunity of the MicroRPM in relation to EMC performance.

The MicroRPM should be used only with the accessories (RS232 cables) supplied (which are referenced in the accessories section of this manual). None of the MicroRPM cables should be extended in length by the user.

If any cables are extended by the user or non-approved accessories are used, this may result in an increased level of emissions or decreased level of immunity in relation to the MicroRPM EMC. None of the MicroRPM accessories should be used with other devices, as this may result in an increased level of emissions or decreased level of immunity in relation to the other devices' EMC.

The MicroRPM has a minimum basic performance – The respiratory pressure readings on the product must remain within a tolerance of +/- 3%, and the unit firmware must not cease responding. Warning: In the event the product is operated in the presence of significant electromagnetic fields (particularly in the frequency range 40-60MHz) while in the PC connected mode, ensure the results on the unit and PC are the same. If the results differ, then relocating the product away from sources of interference should resolve any issue.

WARNING: The MicroRPM should not be used adjacent to or stacked with other equipment. If adjacent or stacked use with other equipment is necessary, the MicroRPM and the other equipment should be observed / monitored, to verify normal operation in the configuration in which it will be used.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions			
The MicroRPM is intended for use in the electromagnetic environment specified below. The customer or the user of the MicroRPM should assure that it is used in such an environment.			
Emission Test	Compliance	Electromagnetic Environment - Guidance	
RF emissions CISPR11	Group 1	The MicroRPM uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF emissions CISPR11	Class B	The MicroRPM is suitable for use in all establishments, including residential establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
Harmonic Emissions IEC 61000-3-2	[Not Applicable]		
Voltage fluctuations / flicker emissions IEC 61000-3-3	[Not applicable]		

Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
The MicroRPM is intended for use in the electromagnetic environment specified below. The customer or the user of the MicroRPM should assure that it is used in such an environment.			
Immunity Test	IEC 60601 test level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6 kV contact +/- 8 kV air	+/- 6 kV contact +/- 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	+/- 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.

Surge IEC 61000-4-5	+/- 1 kV line(s) to line (s) +/- 2 kV line(s) to earth	Not Applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% UT (>95% dip in UT) for 0.5% cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (>95% dip in UT) for 5s	Not Applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the micro RPM requires continued operation during power main interruptions, it is recommended that the micro RPM be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE UT is the a.c. mains voltage prior to application of the test level.			

Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
The MicroRPM is intended for use in the electromagnetic environment specified below. The customer or the user of the MicroRPM should assure that it is used in such an environment.			
Immunity Test	IEC 60601 test level	Compliance Level	Electromagnetic Environment -Guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the micro RPM, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

<p>Conducted RF IEC61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150kHz to 80MHz</p> <p>3 V/m 80MHz to 2.5GHz</p>	<p>3 V</p> <p>3 V/m</p>	<p>Recommended separation distance</p> <p>$d = 1.2 \times \sqrt{P}$</p> <p>$d = 1.2 \times \sqrt{P}$ 80MHz to 800MHz</p> <p>$d = 2.3 \times \sqrt{P}$ 800MHz to 2.5GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from a fixed RF transmitters as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radio, amateur radio, AM and FM radio broadcast and TV broadcast, cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the MicroRPM is used exceeds the applicable RF compliance level above, the MicroRPM should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orientating or relocating the MicroRPM.</p>			
<p>b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m</p>			

Recommended separation distances between portable and mobile RF communications equipment and the MicroRPM

The MicroRPM is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the MicroRPM can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the MicroRPM as recommended below, according to the maximum output power of the communications equipment

Rated Maximum Output Power of Transmitter in Watts (W)	Separation Distance in Meters (m) according to Frequency of Transmitter		
	150 KHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1.2 \times \sqrt{P}$	$d = 1.2 \times \sqrt{P}$	$d = 2.4 \times \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Symbols


Type B device

Consult the instructions for use

Caution: Consult the accompanying documents

Date of Manufacture

Manufacturer

Serial number

Single patient use

Federal U.S. law restricts this device to sale by or on the order of a physician. (Rx only)

Nationally Recognized Test Laboratory (NRTL) officially recognized by the American Occupational Safety and Health Administration (OSHA) for OSHA Electrical safety and compliance.

Specifications

Measurements	Maximum Expiratory Pressure (MEP) Maximum Inspiratory Pressure (MIP) Sniff Nasal Inspiratory Pressure (SNIP)
With Puma™	Maximum Rate of Pressure Development (MRPD) Maximum Rate of Relaxation (MRR) Tau
Operating Pressure:	+/- 300cmH2O (±5PSID)
Burst Pressure:	+/- 700cmH2O (±20PSID)
Accuracy:	+/-3%
Resolution:	1cmH2O
Power Supply:	Single Alkaline 9V 6LR61
Dimensions:	170 x 60 x 26mm
Weight (unit with battery):	175g
Weight (in carry case) :	750g
Operating Temperature :	0 - 40 °C
Operating Humidity:	30% - 90% RH
Storage & Transport temperature:	-20 - +70 ° C
Storage & Transport humidity:	10% - 90% RH

Consumables / Supporting Products

Description

Expiratory Pressure Valve Assembly
Inspiratory Pressure Valve Assembly
Bacterial Filters (50 per box)
Extra Small Nasal Probes (10 per box)
Small Nasal Probes (10 per box)
Medium Nasal Probes (10 per box)
Large Nasal Probes (10 per box)
Nasal Probe Adapter
Calibration Screwdriver
Alkaline PP3 Battery
PeraSafe Sterilising Powder – 81g
(Makes 5 litres of solution)
Pressure Calibration Kit
Interface Cable for PUMA Software

To place an order for consumables/supporting products or for general enquiries please contact CareFusion or contact your local CareFusion dealer.

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VI Customer contact information

UK Customers only

For all Sales Order processing for products, training and spare parts, Service and Technical Support enquiries, please contact the following:

CareFusion UK 232 Ltd
UK Customer Service & Support
The Crescent
Jays Close
Basingstoke
RG22 4BS

Customer Service Sales Enquiries:

Telephone: 01256 388550
Email: micro.uksales@carefusion.com

Factory Repair and Administration Enquiries:

Telephone: 01256 388552
Email: micro.ukservice@carefusion.com

Technical Support Enquiries:

Telephone: 01256 388551
Email: support.rt.eu@carefusion.com

International customers only

For all Sales Order processing for products and Spare parts, Service and Technical Support enquiries, please contact the following:

Carefusion Germany 234 GmbH
Customer Service & Support International
Leibnizstrasse 7
97204 Hoechberg
Germany

Customer Service Sales Enquiries:

Telephone: 0049 931 4972 670

Email: micro.internationalsales@carefusion.com

Factory Repair and Administration Enquiries:

Telephone: 0049 931 4972 867

Email: support.admin.eu@carefusion.com

Technical Support Enquiries:

Telephone: 0049 931 4972 127

Email: support.rt.eu@carefusion.com

Notes:

Notes:

CareFusion UK 232 Ltd.,
The Crescent
Jays Close
Basingstoke RG22 4BS
U.K.

CE
0086

JDE no. 36-MAN1294
Drg no. 064-45 WW (Part 1)
Issue 1.2
January 2014
© CareFusion 2014

carefusion.com/micromedical

CareFusion