CareFusion

Spirometry PC Software

Getting Started Guide

Getting Started Guide – English US

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Contents

I.	Welcome – English	8
A.	System Requirements	8
1.	PC Requirements	8
2.	Operating System Requirements.....	9
B.	Installing Spirometry PC Software	9
C.	Spirometry PC Software Settings.....	10
D.	Hardware Key	11
E.	Creating a Database	12
F.	Support.....	12
G.	Connectivity to GDT.....	13
II.	Welcome – CareFusion USA	17
A.	System Requirements	17
1.	PC Requirements	17
2.	Operating System Requirements.....	17
B.	Installing Spirometry PC Software	18
C.	Spirometry PC Software Settings.....	19
D.	Hardware Key	20
E.	Creating a Database	20
F.	Support.....	20
G.	Connectivity to GDT.....	21
III.	Bienvenue – Français	24
A.	Configuration requise.....	24
1.	Configuration du PC	24
2.	Système d'exploitation	25
B.	Installation du logiciel Spirometry PC Software	25
C.	Paramètres du logiciel Spirometry PC Software	26
D.	Clé électronique	27

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Getting Started Guide

I. Welcome – English

Thank you for choosing Spirometry PC Software. This getting started guide provides step-by-step instructions to install and perform initial setup of the software. Once successfully installed, more comprehensive help is available through the Help menu when using the application or alternatively the PDF manual located on the Spirometry PC Software installation CD.

The system is intended for use with paediatric (4 to 17 years) and adult (18 to 99 years) patients in hospitals, physician offices, laboratories and occupational health testing environments.

A. System Requirements

Spirometry PC software requires certain hardware and software components in order to run properly.

1. PC Requirements

An IBM-compatible PC is required, with hardware that meets or exceeds the following minimum requirements.

- Processor: 800 MHz or above

- RAM: 256 MB
- Free Disk Space: 50 MB
- Video: 800x600, 256 colours. It is recommended that a resolution of at least 1280x1024 be used to enable the full benefits of the multi-window interface.
- At least one free USB port. (An additional USB port will be required for USB enabled spiroometers)

2. Operating System Requirements

Spirometry PC Software will run on the following operating systems:

- Windows XP
- Windows Vista 32bit
- Windows Vista 64bit
- Windows 7 32bit
- Windows 7 64bit

B. Installing Spirometry PC Software

Before you begin, please ensure your computer meets the minimum system requirements and the user installing the software has administrative user rights.

1. Close any other applications that are running.
2. Insert the installation CD into your CD-Rom drive.
3. The setup program should launch automatically, displaying the welcome dialog box, click 'Next' to continue. If the setup program does not launch automatically use Windows explorer to manually select the CD drive then open the file named Spirometry PC Software.exe.
4. The license agreement will be displayed. Please read the complete document and ensure you understand fully before accepting the terms of the license. Click 'Next' to continue.

5. The select destination location dialog will be displayed showing the directory where the Spirometry PC Software will be installed. The default location is
C:\Program Files\CareFusion\Spirometry PC Software
To change this location 'select 'Browse'.
6. Click 'Next' to display the Select program manager group dialog. This will be the location where the Spirometry PC Software launch icon will be placed under the Start Menu. Either a new location can be entered or an existing location used.
7. Click 'Next' to display the installation options dialog containing the following option:
 - Add Spirometry PC Software icon to the desktop.
8. Click 'Next' to start copying the files to your system.
9. Once the files have finished copying to your system you will be advised that the Spirometry PC Software has been successfully installed and that it would be advisable to restart your PC before using the application. Click 'Finish' to complete the installation process.

C. Spirometry PC Software Settings

The 'Settings' dialog will be displayed the first time you enter the Spirometry PC Software. Settings include:

- Regional Settings
This settings option allows you to select the region and language settings.
- General Settings
This settings option allows you to select general settings such as height units, weight units, date format and date separator.
- Print Settings
This settings option allows you to select the settings for the printed reports.

- Spirometry Settings
This settings option allows you to change any spirometry settings.
- Spirometry Indices Settings
These options allow you to select which spirometry indices will be reported within various Spirometry PC Software windows and reports.

D. Hardware Key

After installation the below dialog box is displayed.

Prior to using the Spirometry PC Software, you should insert the Spirometry PC Software hardware key into a free USB port on your PC. This will allow access to the Spirometry PC Software program. The hardware key provides freedom of use allowing the Spirometry PC Software to be installed onto numerous computers but can only be activated using the hardware key.

Please ensure that you keep the hardware key safe, if lost a replacement will need to be purchased.

E. Creating a Database

Clicking on 'File' will enable you to:

- Create your own database by selecting 'File' and then 'New'
- Open an existing database by selecting 'Open' and the database name

F. Support

If you experience any problems while using the Spirometry PC Software, please contact our support team who will be happy to help you.

If linking to either EMIS LV, EMIS PCS, EMIS WEB or SystmOne please use the links below for your connection guide.

EMIS: <http://www.micromedical.co.uk/update/spcs/emis/spcs-emis-installer.exe>

SystmOne: <http://www.micromedical.co.uk/update/spcs/S1/spcs-s1-installer.exe>

G. Connectivity to GDT

SPCS provides a GDT interface for interoperability with a practice electronic medical record system (PAS).


The PAS must be installed and configured to launch SPCS in the SPCS installation folder using the /GDT option.

Example:

If the default settings are used SPCS shall be started as follows:

c:\Program Files\CareFusion\Spirometry PC Software\Spirometry PC Software.exe **/GDT2 (SPCS 2.2 and above)**.

When SPCS is first started after installation from your PAS the below dialogue box will be displayed.

The configuration for Incoming and Outgoing directories can be defined and must match the values set in the PAS.

Messages to GDT PMS

In this section the content of the GDT 8315 and 8316 segments are defined. (8315 refers to the SPCS system, 8316 refers to the PAS system).

Message File Naming

The “Short name for this system” is used to identify files intended for the software to read. Please refer to the GDT specification for further details on file naming conventions.

Attention: the GDT filename **must** start with the “short name”.

Result table text length

Please define the width of result table according to your requirements. The width is given in characters and does not include characters of the GDT header.


Messages to GDT PMS (version 2.1 only)

If the SPCS system is configured for GDT version 2.1 additional options are available.


9206 Character set	Select the appropriate character set.
Send link to report	Together with the Outgoing GDT message a text report will be returned to the PAS. You may embed a link (path) to a WORD or PDF report in the GDT message. The default directory used to transmit this file to the PAS is the Windows Temp directory.
6304 File Content	Use this segment to include a fixed textual definition of the GDT message's content according to your requirements.
Report Filename	Use the drop down menu to select how the report is to be labelled and filed, either by Unique File name, or the file name using the Patients ID and exam date/time. There is also the option to select whether the report is saved as a PDF or Word document.

NOTE:

Whenever started from the PAS, SPCS expects an incoming GDT message. If no such message is created the following error message is displayed:

GDT settings can also be reconfigured after initial installation and configuration from the main tool bar in SPCS by selecting Data and Export Options

After initial installation, to perform a spirometry test on the patient, the patient must be selected from the PAS, once the patient is selected, to perform spirometry select Lung Function from the Practice software and SPCS will be launched for a test to be performed on the patient.

The standard work flow is to perform relaxed spirometry followed by forced spirometry. However if relaxed spirometry is not required, the Forced spirometry option can be selected for immediate testing.

After spirometry has been performed, the results should be saved in SPCS, and SPCS should be closed. If the user does not Save and Exit SPCS correctly, results will be exported to the Patient records. When SPCS is closed the spirometry results will automatically be transferred to the patient's records in the Practice system, along with any changes made to the patient demographics.

The indices (parameters) exported will be those that had previously been selected in SPCS using the Customise option and Spirometry Indices settings.

CareFusion

Spirometry PC Software

Getting Started Guide

Federal (USA) law restricts this device to sale by or on the order of a physician or licensed practitioner.

CareFusion,
22745 Savi Ranch Parkway,
Yorba Linda
CA 92887-4668
USA

II. Welcome – CareFusion USA

Thank you for choosing Spirometry PC Software. This getting started guide provides step-by-step instructions to install and perform initial setup of the software. Once successfully installed, more comprehensive help is available through the Help menu when using the application or alternatively the PDF manual located on the Spirometry PC Software installation CD.

A. System Requirements

Spirometry PC software requires certain hardware and software components in order to run properly.

1. PC Requirements

An IBM-compatible PC is required, with hardware that meets or exceeds the following minimum requirements.

- Processor: 800 MHz or above
- RAM: 256 MB
- Free Disk Space: 50 MB
- Video: 800x600, 256 colors. It is recommended that a resolution of at least 1280x1024 be used to enable the full benefits of the multi-window interface.
- At least one free USB port. (An additional USB port will be required for USB enabled spirometers)

2. Operating System Requirements

Spirometry PC Software will run on the following operating systems:

- Windows XP
- Windows Vista 32bit
- Windows Vista 64bit
- Windows 7 32bit
- Windows 7 64bit

B. Installing Spirometry PC Software

Before you begin, please ensure your computer meets the minimum system requirements and the user installing the software has administrative user rights.

1. Close any other applications that are running.
2. Insert the installation CD into your CD-Rom drive.
3. The setup program should launch automatically, displaying the welcome dialog box, click 'Next' to continue. If the setup program does not launch automatically use Windows explorer to manually select the CD drive then open the file named Spirometry PC Software.exe.
4. The license agreement will be displayed. Please read the complete document and ensure you understand fully before accepting the terms of the license. Click 'Next' to continue.
5. The select destination location dialog will be displayed showing the directory where the Spirometry PC Software will be installed. The default location is
C:\Program Files\Cardinal Health\Spirometry PC Software
To change this location 'select 'Browse'.
6. Click 'Next' to display the Select program manager group dialog. This will be the location where the Spirometry PC Software launch icon will be placed under the Start Menu. Either a new location can be entered or an existing location used.
7. Click 'Next' to display the installation options dialog containing the following option:
 - Add Spirometry PC Software icon to the desktop.
8. Click 'Next' to start copying the files to your system.

- Once the files have finished copying to your system you will be advised that the Spirometry PC Software has been successfully installed and that it would be advisable to restart your PC before using the application. Click 'Finish' to complete the installation process.

C. Spirometry PC Software Settings

The 'Settings' dialog will be displayed the first time you enter the Spirometry PC Software. Settings include:

- Regional Settings

This settings option allows you to select the region and language settings.

- General Settings

This settings option allows you to select general settings such as height units, weight units, date format and date separator.

- Print Settings

This settings option allows you to select the settings for the printed reports.

- Spirometry Settings


This settings option allows you to change any spirometry settings.

- Spirometry Indices Settings

These options allow you to select which spirometry indices will be reported within various Spirometry PC Software windows and reports.

D. Hardware Key

After installation the below dialog box is displayed.

Prior to using the Spirometry PC Software, you should insert the Spirometry PC Software hardware key into a free USB port on your PC. This will allow access to the Spirometry PC Software program. The hardware key provides freedom of use allowing the Spirometry PC Software to be installed onto numerous computers but can only be activated using the hardware key.

Please ensure that you keep the hardware key safe, if lost a replacement will need to be purchased.

E. Creating a Database

Clicking on 'File' will enable you to:

- Create your own database by selecting 'File' and then 'New'
- Open an existing database by selecting 'Open' and the database name

F. Support

If you experience any problems while using the Spirometry PC Software, please contact our support team who will be happy to help you.

G. Connectivity to GDT

SPCS provides a GDT interface for interoperability with a practice electronic medical record system (PAS).


The PAS must be installed and configured to launch SPCS in the SPCS installation folder using the /GDT option.

Example:

If the default settings are used SPCS shall be started as follows:

c:\Program Files\CareFusion\Spirometry PC Software\Spirometry PC Software.exe **/ GDT2 (SPCS 2.2 and above)**.

When SPCS is first started after installation from your PAS the below dialogue box will be displayed.

The configuration for Incoming and Outgoing directories can be defined and must match the values set in the PAS.

Messages to GDT PMS

In this section the content of the GDT 8315 and 8316 segments are defined. (8315 refers to the SPCS system, 8316 refers to the PAS system).

Message File Naming

The "Short name for this system" is used to identify files intended for the software to read. Please refer to the GDT specification for further details on file naming conventions.

Attention: the GDT filename **must** start with the "short name".

Result table text length

Please define the width of result table according to your requirements. The width is given in characters and does not include characters of the GDT header.


Messages to GDT PMS (version 2.1 only)

If the SPCS system is configured for GDT version 2.1 additional options are available.


9206 Character set	Select the appropriate character set.
Send link to report	Together with the Outgoing GDT message a text report will be returned to the PAS. You may embed a link (path) to a WORD or PDF report in the GDT message. The default directory used to transmit this file to the PAS is the Windows Temp directory.
6304 File Content	Use this segment to include a fixed textual definition of the GDT message's content according to your requirements.
Report Filename	Use the drop down menu to select how the report is to be labelled and filed, either by Unique File name, or the file name using the Patients ID and exam date/time. There is also the option to select whether the report is saved as a PDF or Word document.

NOTE:

Whenever started from the PAS, SPCS expects an incoming GDT message. If no such message is created the following error message is displayed:

GDT settings can also be reconfigured after initial installation and configuration from the main tool bar in SPCS by selecting Data and Export Options

After initial installation, to perform a spirometry test on the patient, the patient must be selected from the PAS, once the patient is selected, to perform spirometry select Lung Function from the Practice software and SPCS will be launched for a test to be performed on the patient.

The standard work flow is to perform relaxed spirometry followed by forced spirometry. However if relaxed spirometry is not required, the Forced spirometry option can be selected for immediate testing.

After spirometry has been performed, the results should be saved in SPCS, and SPCS should be closed. If the user does not Save and Exit SPCS correctly, results will be exported to the Patient records. When SPCS is closed the spirometry results will automatically be transferred to the patient's records in the Practice system, along with any changes made to the patient demographics.

The indices (parameters) exported will be those that had previously been selected in SPCS using the Customise option and Spirometry Indices settings.

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Customer contact information

UK Customers only

For all Sales Order processing for products, training and spare parts, Service and Technical Support enquiries, please contact the following:

CareFusion UK 232 Ltd
UK Customer Service & Support
The Crescent
Jays Close
Basingstoke
RG22 4BS

Customer Service Sales Enquiries:

Telephone: 01256 388550
Email: micro.uksales@carefusion.com

Factory Repair and Administration Enquiries:

Telephone: 01256 388552
Email: micro.ukservice@carefusion.com

Technical Support Enquiries:

Telephone: 01256 388551
Email: support.rt.eu@carefusion.com

International customers only

For all Sales Order processing for products and Spare parts, Service and Technical Support enquiries, please contact the following:

Carefusion Germany 234 GmbH
Customer Service & Support International
Leibnizstrasse 7
97204 Hoechberg
Germany

Customer Service Sales Enquiries:

Telephone: 0049 931 4972 670

Email: micro.internationalsales@carefusion.com

Factory Repair and Administration Enquiries:

Telephone: 0049 931 4972 867

Email: support.admin.eu@carefusion.com

Technical Support Enquiries:

Telephone: 0049 931 4972 127

Email: support.rt.eu@carefusion.com

CareFusion UK 232 Ltd.,
The Crescent
Jays Close
Basingstoke
RG22 4BS

CE
0086

JDE no. 36-MAN1284
Drg no. 107-162WW
Issue 1.4
October 2013
© CareFusion 2013

carefusion.com/micromedical

CareFusion